



---

## Contents of Kit:

### 48 Reactions

Component	Size	Number
Reaction Pellets	48 pellets	6 bottles
Pre-filled 100 µL PD1 Pellet diluent tubes	48 tubes	6 strips
AMP1 extraction buffer	55 mL	1 bottle
Sample extraction bags, mesh	50 bags	1 pack

### Questions or Technical Support:

Phone: 800-622-4342 (toll-free) or 574-264-2014

Fax: 574-264-2153

E-mail: [info@agdia.com](mailto:info@agdia.com) for sales and general product information  
[techsupport@agdia.com](mailto:techsupport@agdia.com) for technical information and troubleshooting

Web: [www.agdia.com](http://www.agdia.com)

AmplifyRP Test Kits employ recombinase polymerase amplification (RPA) technology, developed by TwistDx Limited, U.K. Use of the RPA process and probe technologies are protected by US patents 7,270,981 B2, 7,399,590 B2, 7,435,561 B2, 7,485,428 B2 and foreign equivalents in addition to pending patents.

AmplifyRP® is a registered trademarks of Agdia, Inc.