

Tech Note: Sample Extraction Buffer Change

SRA 36600 • *Angelonia flower break virus* (AnFBV)

Technical Note

Based on customer feedback, Agdia has been pursuing an initiative to simplify the extraction process as much as possible across our ELISA product line. The extraction buffer for this test will be changing from General Extract Buffer 3 (GEB3) to General Extract Buffer (GEB) starting March 1st, 2021. By making this change, we seek to eliminate the need for duplicate sampling and extraction events when using more than one ELISA test with each sample.

Please note that GEB extraction buffer must now be purchased in addition to the reagent set. Since GEB is not a proprietary buffer formulation, you may also opt to make it from raw chemicals sourced outside of Agdia. While GEB3 will no longer be provided with the test listed above, it will remain available for purchase as an optional accessory and use of it remains fully supported by Agdia.

Sample Extraction Buffer Change (Beginning March 1 st , 2021)		
	Previous Protocol	New Protocol
Sample Extraction Buffer	GEB3 (ACC 00360)	GEB (ACC 00955)
Included with reagent set?	Yes	No (Purchased separately)

For answers to your technical questions, please contact us at techsupport@agdia.com.



Agdia, Inc.
52642 County Road 1
Elkhart, IN 46514
574-264-2014 / 800-622-4342
www.agdia.com / info@agdia.com

n221.1
Revised: 02/25/2021
Page 1 of 1