

Tech Note: Product Change

PSA/SRA 20501 • *Impatiens necrotic spot virus* (INSV) ELISA

Technical Note

We have recently updated the antibodies for *Impatiens necrotic spot virus* (INSV). These updates have resulted in the following changes outlined below. Please visit the product webpage for the [validation report](#). Previous antibodies and components are not interchangeable with the new version.

Changeover Guide

Do not use previous test components with the new test components.

| Item | Previous Assay | New Assay |
|-----------------------------|-----------------|---------------------------|
| Test Format | Compound-ELISA | DAS-ELISA |
| PathoScreen® Kit User Guide | m20.7 | m14.8 |
| Reagent Set User Guide | m12.5 | m16.6 |
| Capture Antibody | 00493 and below | 00494 and above |
| Detection Antibody | 00540 and below | Not used in current assay |
| Enzyme Conjugate | 00540 and below | 00541 and above |
| Coated Plates | 00005 and below | 00006 and above |
| Validation Report Number | p284.1 | p345 |

For answers to your technical questions, please contact us at techsupport@agdia.com.

PathoScreen® is a registered trademark of Agdia, Inc.



Agdia, Inc.
52642 County Road 1
Elkhart, IN 46514
574-264-2014 / 800-622-4342
www.agdia.com / info@agdia.com

n243
Revised: 11/20/2023
Page 1 of 1