


# Tech Note: Product Change

PSA/SRA 36600 • *Angelonia flower break virus* (AnFBV) ELISA



## Technical Note

We have recently updated the antibodies for *Angelonia flower break virus* (AnFBV). These updates have resulted in the following changes outlined below. Please visit the product webpage for the [validation report](#). Previous antibodies and components are not interchangeable with the new version.

## Changeover Guide

Do not use previous test components with the new test components.

Item	Previous Assay	New Assay
Test Format	DAS-ELISA	DAS-ELISA
PathoScreen® Kit User Guide	N/A	m404.1
Reagent Set User Guide	m392.3	m392.3
Capture Antibody	00300 and below	00301 and above
Enzyme Conjugate	00351 and below	00352 and above
Coated Plates	N/A	00011 and above
Validation Report Number	p147.2	p351

For answers to your technical questions, please contact us at [techsupport@agdia.com](mailto:techsupport@agdia.com).

PathoScreen® is a registered trademark of Agdia, Inc.



Agdia, Inc.  
52642 County Road 1  
Elkhart, IN 46514  
574-264-2014 / 800-622-4342  
[www.agdia.com](http://www.agdia.com) / [info@agdia.com](mailto:info@agdia.com)

n245  
Revised: 04/18/2024  
Page 1 of 1