


Tech Note: Product Change

PSA/SRA 31505 • *Plum pox virus* (PPV) ELISA



Technical Note

We have recently updated the antibodies for *Plum pox virus* (PPV) and now offer a PathoScreen® Kit with a precoated plate. These updates have resulted in the following changes outlined below. Please visit the product webpage for the [validation report](#). Previous antibodies and components are not interchangeable with the new version.

Changeover Guide

Do not use previous test components with the new test components.

| Item | Previous Assay | New Assay |
|-----------------------------|-----------------|-----------------|
| Test Format | DAS-ELISA | DAS-ELISA |
| PathoScreen® Kit User Guide | N/A | m471 |
| Reagent Set User Guide | m243.4 | m243.5 |
| Capture Antibody | 01039 and below | 01040 and above |
| Enzyme Conjugate | 01063 and below | 01064 and above |
| Coated Plates | N/A | 00028 and above |
| Validation Report Number | p117.2 | p357 |

For answers to your technical questions, please contact us at techsupport@agdia.com.

PathoScreen® is a registered trademark of Agdia, Inc.



Agdia, Inc.
52642 County Road 1
Elkhart, IN 46514
574-264-2014 / 800-622-4342
www.agdia.com / info@agdia.com

n247
Revised: 10/08/2024
Page 1 of 1