



Technical Note and Alert

The ImmunoStrip for Citrus tristeza virus (CTV) has been redeveloped to improve product performance. The new version makes it easier for users to interpret low-level positives, resulting in lower risk of non-specific results. Validation studies indicate that specificity and sensitivity characteristics are comparable between both versions.

Changeover Guide

Do not use previous test components with the new test components.

Item	Previous ImmunoStrip	New ImmunoStrip
ImmunoStrips	00026 and below	00027 and above
User Guide	m264.1	m264.1
Validation Report	p108.1	p108.2

For answers to your technical questions, please contact us at techsupport@agdia.com.