


Tech Note: Product Change

SRA 88001 • *Blueberry leaf mottle virus* (BLMoV) ELISA



Technical Note

We have recently updated the antibodies for *Blueberry leaf mottle virus* (BLMoV). The test system has changed from a three bottle triple antibody sandwich system (TAS-ELISA) to a two bottle double antibody sandwich system (DAS-ELISA). See the [user guide](#), m329.3, for protocol details. These updates have resulted in the following changes outlined below. Please visit the product webpage for the [validation report](#). Previous antibodies and components are not interchangeable with the new version.

Changeover Guide

Do not use previous test components with the new test components.

Item	Previous Assay	New Assay
Test Format	TAS-ELISA	DAS-ELISA
Reagent Set User Guide	m12.5	m329.3
Capture Antibody	00471 and below	00472 and above
Detection Antibody	00469 and below	N/A
Enzyme Conjugate	00469 and below	00470 and above
Validation Report Number	p204.1	p204.2

For answers to your technical questions, please contact us at techsupport@agdia.com.



Agdia, Inc.
52642 County Road 1
Elkhart, IN 46514
574-264-2014 / 800-622-4342
www.agdia.com / info@agdia.com

n252
Revised: 02/25/2025
Page 1 of 1