



Technical Note and Alert

The ImmunoStrip® for *Melon severe mosaic virus* (MeSMV) has been redeveloped to ensure an uninterrupted supply of this product for years to come. The new version has superior performance characteristics, including the ability to detect lower viral concentrations compared to the previous version. The extraction buffer has changed from SEB9 (previous version) to SEB1 (new version), which is compatible with the majority of Agdia's pathogen ImmunoStrips®.

Changeover Guide

Do not use previous test components with the new test components.

Item	Previous ImmunoStrip	New ImmunoStrip
ImmunoStrips	00006 and below	00007 and above
Extraction Buffer	SEB9	SEB1
User Guide	m325.2	m175.6
Validation Report	p228.3	p368

For answers to your technical questions, please contact us at techsupport@agdia.com.