



Technical Note and Alert

The extraction buffer for the Melon severe mosaic virus (MeSMV) ImmunoStrip® has been changed to ensure optimal performance of this product. The extraction buffer has changed from SEB1 (previous version) to SEB4 (new version). The SEB4 extraction buffer has significantly reduced false positive reactions. No other ImmunoStrip® components have been changed. The diagnostic and analytical sensitivity have not been affected. If you have unused lot 00007 ImmunoStrips that were received with SEB1 extraction buffer, please reach out to receive a replacement of SEB4 extraction buffer.

Changeover Guide

Do not use previous test components with the new test components.

Item	ImmunoStrip before 4/2/2026	ImmunoStrip after 4/2/2026
ImmunoStrips	00007	00007 and above
Extraction Buffer	SEB1	SEB4
User Guide	m175.6	m494
Validation Report	p368	p368.1

For answers to your technical questions, please contact us at techsupport@agdia.com.