


Content List: DAS-ELISA Reagent Sets

Proprietary Buffers • Components and Sizes



96 Testwells

Component	Size	Number
Capture antibody	0.150 mL	1 bottle
Peroxidase enzyme conjugate	0.150 mL	1 bottle
TMB peroxidase substrate solution, 1X	25 mL	1 bottle
MRS component, 5X	3 mL	1 bottle
96-well microtiter plate, strip	N/A	1 plate

500 Testwells

Component	Size	Number
Capture antibody	0.275 mL	1 bottle
Peroxidase enzyme conjugate	0.275 mL	1 bottle
TMB peroxidase substrate solution, 1X	60 mL	1 bottle
MRS component, 5X	11 mL	1 bottle
96-well microtiter plate, solid	N/A	5 plates

1000 Testwells

Component	Size	Number
Capture antibody	0.525 mL	1 bottle
Peroxidase enzyme conjugate	0.525 mL	1 bottle
TMB peroxidase substrate solution, 1X	60 mL	2 bottles
MRS component, 5X	21 mL	1 bottle
96-well microtiter plate, solid	N/A	10 plates

5000 Testwells

Component	Size	Number
Capture antibody	2.525 mL	1 bottle
Peroxidase enzyme conjugate	2.525 mL	1 bottle
TMB peroxidase substrate solution, 1X	550 mL	1 bottle
MRS component, 5X	102 mL	1 bottle
96-well microtiter plate, solid	N/A	50 plates

For answers to your technical questions, please contact us at techsupport@agdia.com.



Agdia, Inc.
52642 County Road 1
Elkhart, IN 46514
574-264-2014 / 800-622-4342
www.agdia.com / info@agdia.com

CL_1_SRP
Revised: 10/25/2021
Page 1 of 1