



---

## Contents of Kit:

### 48 Reactions

Component	Size	Number
Reaction Pellets	48 pellets	6 bottles
Pre-filled 100 µL PD2 Pellet diluent tubes	48 tubes	6 strips
AMP1 extraction buffer	55 mL	1 bottle
Microcentrifuge tubes	1.5 mL	50 tubes

### Questions or Technical Support:

Phone: 800-622-4342 (toll-free) or 574-264-2014

Fax: 574-264-2153

E-mail: [info@agdia.com](mailto:info@agdia.com) for sales and general product information

[techsupport@agdia.com](mailto:techsupport@agdia.com) for technical information and troubleshooting

Web: [www.agdia.com](http://www.agdia.com)

AmplifyRP® is a registered trademarks of Agdia, Inc.