



Contents of Kit:

48 Reactions

Component	Size	Number
Reaction Pellets	48 pellets	6 bottles
Pre-filled 100 µL PD1 Pellet diluent tubes	48 tubes	6 strips
GEB2-filled mesh extraction bags	50 bags	2 packs
PCR tubes, 0.2 mL	50 tubes	2 packs

Questions or Technical Support:

Phone: 800-622-4342 (toll-free) or 574-264-2014

Fax: 574-264-2153

E-mail: info@agdia.com for sales and general product information
techsupport@agdia.com for technical information and troubleshooting

Web: www.agdia.com

AmplifyRP® is a registered trademarks of Agdia, Inc.